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





Precision Technical Specifications

- Cartridge Mixing technique (XL, L): Use a 1:1/1:2 dispensing cartridge gun, with a yellow static mixing tip (4.2-12D Diameter).
- Cartridge Mixing technique (M & H): Use a 1:1/1:2 dispensing cartridge gun, with a blue static mixing tip (6.5-11D Diameter), or 5:1 automix unit and green dynamic mixing tip for FatPack™ foil bags.

Precision Extra Lite Body Color: Orange

- Linear dimensional change: max. -0.1 %
- Recovery from deformation: min. 99.6 %
- Strain in compression: <4.5 %
- Shore A Hardness: 44-49

Precision Lite Body Color: Green

- Linear dimensional change: max. -0.1 %
- Recovery from deformation: min. 99.6 %
- Strain in compression: <4.0 %
- Shore A Hardness: 46-51

Precision Medium Body Color: Blue

- Linear dimensional change:
- Recovery from deformation:
- Strain in compression:
- Shore A Hardness:

Cartridge

- max. -0.2 %
- min. 99.6 %
- <3.5 %
- 58-62

FatPak

- max. -0.25%
- min. 99.6%
- <4%
- 56-62

Precision Heavy Body Color: Purple

- Linear dimensional change:
- Recovery from deformation:
- Strain in compression:
- Shore A Hardness:

Cartridge

- max. -0.2%
- min. 99.3 %
- max. 3.5%
- 59-63

FatPak

- max. -0.25%
- min. 99.3%
- max. 3.5%
- 57-63

Precision Putty & Putty Paks Color: Purple

- Linear dimensional change: max. -0.2%
 - Recovery from deformation: min. 99.0%
 - Strain in compression: <2.5%
 - Shore A Hardness: 68-72
-
- Storage conditions: In a dry place, at 18 –25 °C (64 –77 °F) and relative humidity of 50% or less. Do not expose to direct heat or sun. Store cartridges upright.
 - Disinfection: Use a cold sterilant.
 - Electroplating: Impressions may be copper or silver plated.
 - Casting: The impression can be poured after 1 hour.

Indications for Use

Precision is an addition-reaction silicon impression material that is ideal for dual viscosity impression techniques both one-step and two-step as well as for a variety of dental procedures including, but not limited to, crowns, bridges, inlays, partial and complete dentures, denture repairs, implants and for forming functional peripheries.


Working and Setting Times

Precision offers a variable working time and setting reaction that is temperature controlled. If placed in the mouth immediately after extrusion, Precision will allow 45 seconds of intra-oral working time. If more time is required, Precision will allow for 1:45 extra-oral working time. Therefore, depending on how much working time is required (45 seconds – 1:45), the total set time will range from 2:45 – 3:45 minutes.

Directions for Use


Instructions for cartridge:

1. Insert the cartridge into the cartridge dispensing gun.
2. Remove the cartridge closure.
3. Before applying the mixing tip, extrude impression material until it exits both chambers at the same rate.
4. Attach the mixing tip on the cartridge and extrude the required amount of material directly onto the prepared tooth or into an impression tray or delivery syringe.
5. Leave the filled mixing tip on the cartridge after use.
6. Before using the cartridge again, remove and discard the used mixing tip. Before re-applying another mixing tip, be sure to extrude material so the product is exiting both chambers at the same rate.
7. Attach a fresh mixing tip and proceed as usual.
8. The cartridge dispensing gun can be cold sterilized.


Tray Filling Time:

Up to 1:45 min.
(73.4° F / 23° C)

Intraoral Syringe Time:

Up to 1 min.
(95° F / 35° C)

Total Set Time:

2:45 – 3:45 min.
(73.4° F / 23° C)

Instructions for FatPak™

Please refer to the illustrations in the back of this booklet.

1. Before first use, turn the white safety pin on the activating head of the large foil bag in the direction indicated by the arrow in the drawing, then remove and dispose of the pin (fig. 1).
2. Push both foil bags into the FatPak cartridge body. Be sure that the cartridge body and activating head notches are aligned. This will ensure optimum alignment of the outlet nozzles (fig. 2).
3. Press the activating heads onto the cartridge body until they lock firmly into final position (fig. 3 & 4). The foil bags are automatically pierced by a pin when the activating heads are pressed onto the cartridge body.
4. Insert the filled cartridge body into the automix unit.
5. Switch the device on and bleed about 1/4 inch of material to ensure even flow of base and catalyst.
6. Assemble the dynamic mixer tip and close the locking clip for the mixer located on the cartridge body. Note: Use only DenMat® FatPak dynamic mixer tips for optimum performance.

Instructions for Putty Jars

1. Scoop even amounts of base and catalyst from jars.
2. Make sure to keep base and catalyst scoops separate from each other.
3. Knead putty thoroughly until the color is homogenous and streak free.
Note: Do not mix while wearing latex gloves: vinyl or nitrile gloves are recommended.
4. Load mixed putty into tray.

Instructions for Putty Paks®

1. Peel off lid from Putty Pak.
2. Fold base and catalyst together.
3. Squeeze product together.
4. Remove product from Pak and mix thoroughly until color is homogenous and streak free.
Note: Do not mix while wearing latex gloves: vinyl or nitrile gloves are recommended.
5. Load mixed putty into tray.

Special Notes & Precautions

- Cartridge material must be used with a 1:1/1:2 dispensing cartridge gun system. FatPak material can be used in a 5:1 compatible automix unit with dynamic mixer tips.
- Solutions applied prior to impression taking (such as retraction solution) may impede the setting reaction of the impression. Therefore, material must be cleaned off thoroughly with water and dried.
- Do not leave remains of the impression material in the sulcus.
- Avoid splattering this material onto clothing as it cannot be removed by dry cleaning.
- The possibility of allergic reactions to the product in sensitive persons cannot be fully excluded. Possible reactions may include temporary numbness or tingling in the soft tissue.
- When taking an impression adjacent to an existing fixed bridge, make sure all undercuts are sufficiently blocked out.
- All intraoral areas touched by latex gloves should be rinsed and dried prior to impression taking.

The information provided by DenMat is based on comprehensive research and experience in application technology. Results are furnished to the best of our knowledge, subject to technical changes within the framework of product development. However, users must comply with and consider all recommendations and information in connection with any use.

Compatibility

Precision is an addition-reaction type silicone impression material and should not be combined with standard or condensation type silicon materials or polyether materials.

For professional use only. U.S.A. Federal Law restricts this device to sale by or on the order of a dentist.

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